

# Product Datasheet

## Murine IL-2 ELISA Set 861.000

Product Name	Murine IL-2 ELISA Set
Catalog Number	861.000
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	15.6-500 pg/ml
Sensitivity	5 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes both natural and recombinant murine IL-2
Sample Types	Serum, Cell culture supernatant