

Product Datasheet

Murine IL-6 ELISA Set 861.020

Product Name	Murine IL-6 ELISA Set
Catalog Number	861.020
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	15.6-500 pg/ml
Sensitivity	10 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes both natural and recombinant murine IL-6
Sample Types	Serum, Cell culture supernatant