

# Product Datasheet

## Murine IL-10 ELISA Set 861.030

Product Name	Murine IL-10 ELISA Set
Catalog Number	861.030
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	31.25-1000 pg/ml
Sensitivity	20 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes both natural and recombinant murine IL-10
Sample Types	Serum, Cell culture supernatant