

## Product Datasheet

### Murine IFN gamma ELISA Set 861.050

<b>Product Name</b>	Murine IFN gamma ELISA Set
<b>Catalog Number</b>	861.050
<b>Description</b>	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
<b>Range</b>	31.25-1000 pg/ml
<b>Sensitivity</b>	15 pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	IFN-g, IFN-gamma
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other murine cytokines
<b>Specificity</b>	Recognizes both natural and recombinant murine IFN-g
<b>Sample Types</b>	Serum, Cell culture supernatant