

Product Datasheet

Murine IL-2 ELISpot Set 862.001

Product Name	Murine IL-2 ELISpot Set
Catalog Number	862.001
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent. Spots read visually or with reader.
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes natural murine IL-2
References	Nishimura, Y. et al., Immunology, 2002; 107(2): 190-8. - Pubmed link