

Product Datasheet

Murine IL-2 ELISpot Set 862.001.010

Product Name	Murine IL-2 ELISpot Set
Catalog Number	862.001.010
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, and blocking reagent.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes natural murine IL-2