

Product Datasheet

Rat TNF alpha ELISA Kit 865.000.192

Product Name	Rat TNF alpha ELISA Kit
Catalog Number	865.000.192
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	15 pg/ml
Tested Applications	ELISA
Alternative Names	TNFa
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes both natural and recombinant rat TNF-a
Sample Types	Serum, Cell culture supernatant