

Product Datasheet

Rat IFN gamma ELISA Kit 865.010.048

Product Name	Rat IFN gamma ELISA Kit
Catalog Number	865.010.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	10 pg/ml
Tested Applications	ELISA
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes both natural and recombinant rat IFN-g
Sample Types	Serum, Cell culture supernatant