

# Product Datasheet

## Rat IL-4 ELISA Kit 865.020.096

Product Name	Rat IL-4 ELISA Kit
Catalog Number	865.020.096
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	15 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes both natural and recombinant rat IL-4
Sample Types	Serum, Cell culture supernatant