

## Product Datasheet

### Rat IFN gamma ELISpot Set 867.011.015

|                     |                                                                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name        | Rat IFN gamma ELISpot Set                                                                                                                            |
| Catalog Number      | 867.011.015                                                                                                                                          |
| Description         | Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, and blocking reagent. |
| Tested Applications | ELISpot                                                                                                                                              |
| Alternative Names   | IFN-g, IFN-gamma                                                                                                                                     |
| Brand               | Diaclone                                                                                                                                             |
| Storage             | Store at +4°C                                                                                                                                        |
| Note                | Spots read visually or with reader                                                                                                                   |
| Cross Reactivity    | No cross reactivity with other rat cytokines                                                                                                         |
| Specificity         | Recognizes natural rat IFN-g                                                                                                                         |