

Product Datasheet

Human IL-10 ELISpot Pair 869.030.010

Product Name	Human IL-10 ELISpot Pair
Catalog Number	869.030.010
Description	Diaclone ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10.
Specificity	Recognizes natural human IL-10