

# Product Datasheet

## Human TNF alpha ELISpot Pair 869.050.010

Product Name	Human TNF alpha ELISpot Pair
Catalog Number	869.050.010
Description	Diaclone ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
Tested Applications	ELISpot
Alternative Names	TNFa
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with soluble TNFR I and TNFR II.
Specificity	Recognizes natural human TNF-a