

# Product Datasheet

## Human Perforin ELISpot Pair 869.140.010

Product Name	Human Perforin ELISpot Pair
Catalog Number	869.140.010
Description	Diaclone ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes natural human Perforin