

# Product Datasheet

## Human IL-17A/F ELISpot Pair 869.170.010

Product Name	Human IL-17A/F ELISpot Pair
Catalog Number	869.170.010
Description	Diaclone ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes natural human IL-17A/F