

Product Datasheet

Murine IFN gamma ELISpot Pair 870.050.010

Product Name	Murine IFN gamma ELISpot Pair
Catalog Number	870.050.010
Description	Diaclone ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
Tested Applications	ELISpot
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes natural murine IFN-g