

Product Datasheet

Rat IFN gamma ELISpot Pair 871.020.010

Product Name	Rat IFN gamma ELISpot Pair
Catalog Number	871.020.010
Description	Diaclone ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
Tested Applications	ELISpot
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes natural rat IFN-g