

# Product Datasheet

## Rat TNF alpha ELISpot Pair 871.030.010

Product Name	Rat TNF alpha ELISpot Pair
Catalog Number	871.030.010
Description	Diaclone ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
Tested Applications	ELISpot
Alternative Names	TNFa
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes natural rat TNF-a