

# Product Datasheet

## Rat IFN gamma ELISA Set 872.000

Product Name	Rat IFN gamma ELISA Set
Catalog Number	872.000
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	31.25-1000 pg/ml
Sensitivity	10 pg/ml
Tested Applications	ELISA
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes both natural and recombinant rat IFN-g
Sample Types	Serum, Cell culture supernatant