

Product Datasheet

Rat TNF alpha ELISA Set 872.010

Product Name	Rat TNF alpha ELISA Set
Catalog Number	872.010
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	31.25-1000 pg/ml
Sensitivity	15 pg/ml
Tested Applications	ELISA
Alternative Names	TNFa
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes both natural and recombinant rat TNF-a
Sample Types	Serum, Cell culture supernatant