

# Product Datasheet

## Rat IL-4 ELISA Set 872.020

|                     |  |
|---------------------|--|
| Product Name        | Rat IL-4 ELISA Set   |
| Catalog Number      | 872.020  |
| Description         | Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition. |
| Range               | 31.25-1000 pg/ml   |
| Sensitivity         | 15 pg/ml   |
| Tested Applications | ELISA  |
| Brand               | Diaclone   |
| Storage             | Store at +4°C  |
| Cross Reactivity    | No cross reactivity with other rat cytokines   |
| Specificity         | Recognizes both natural and recombinant rat IL-4   |
| Sample Types        | Serum, Cell culture supernatant  |