

# Product Datasheet

## Rat IL-4 ELISA Set 872.020

Product Name	Rat IL-4 ELISA Set
Catalog Number	872.020
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	31.25-1000 pg/ml
Sensitivity	15 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes both natural and recombinant rat IL-4
Sample Types	Serum, Cell culture supernatant