

Product Datasheet

Human IL-15 ELISA Kit 873.000

Product Name	Human IL-15 ELISA Kit
Catalog Number	873.000
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	18.3 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reaction with IL-2, IL-4, IL-7, IL-9, IL-10, IL-13, IL-21 and IFN-g
Specificity	Recognizes both natural and recombinant human IL-15
Sample Types	Serum, Plasma, Cell culture supernatant