

Product Datasheet

Human IL-31 ELISA Kit 873.010

Product Name	Human IL-31 ELISA Kit
Catalog Number	873.010
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	18.3 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with IL-1b, IL-2, IL-4, IL-6, IL-6R, Gp130, IL-17F, IL-27, IL-33, GM-CSF
Specificity	Recognizes both natural and recombinant human IL-31
Sample Types	Serum, Plasma and other body fluids, Cell culture supernatants, Buffered solutions