

Product Datasheet

Human IL-33 ELISA Kit 873.020

Product Name	Human IL-33 ELISA Kit
Catalog Number	873.020
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	12.2 pg/ml
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with IL-1b, IL-2, IL-1a, IL-10, IL-12p70, Gp130, IL-17A, IL-23, IFNg, TNFa
Specificity	Recognizes both natural and recombinant human IL-33
Sample Types	Serum, Plasma and other body fluids, Cell culture supernatants, Buffered solutions