

Product Datasheet

Human GM-CSF ELISA Kit 873.040

Product Name	Human GM-CSF ELISA Kit
Catalog Number	873.040
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	15.6-500 pg/ml
Sensitivity	4.8 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant human GM-CSF
Sample Types	Serum, Plasma, Cell culture supernatant