

# Product Datasheet

## Human VEGF165 ELISA Kit 873.070

Product Name	Human VEGF165 ELISA Kit
Catalog Number	873.070
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	15.9pg/ml
Tested Applications	ELISA
Alternative Names	VEGF-A
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reaction with VEGF-121, VEGF-145, VEGF-B, VEGF-C, VEGF-D, VEGF-R1, Ang-2. Shows cross reactivity with VEGF-189 isoform
Specificity	Recognizes human VEGF165 isoforms
Sample Types	Serum, Plasma, Cell culture supernatant