

Product Datasheet

Human Angiopoietin-2 ELISA Kit 873.080

Product Name	Human Angiopoietin-2 ELISA Kit
Catalog Number	873.080
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	125-4000 pg/ml
Sensitivity	25 pg/ml
Tested Applications	ELISA
Alternative Names	ANG-2,ANGPT-2
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with VEGF-A, GM-CSF, IL-34, IL-29 IFNg, IL-2, IL-10, CD138, MICB, MCP-1
Specificity	Recognizes natural and recombinant Angiopoietin-2
Sample Types	Serum, Cell culture supernatant