

# Product Datasheet

## Human IFN $\gamma$ / IL-10 Dual Fluorospot Set 874.012

|                     |                                                                                                                                                                                                                                                                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name        | Human IFN $\gamma$ / IL-10 Dual Fluorospot Set                                                                                                                                                                                                                                                                                                 |
| Catalog Number      | 874.012                                                                                                                                                                                                                                                                                                                                        |
| Description         | Diaclone Dual Fluorospot Sets include capture antibody for cytokine 1, FITC-conjugated detection antibody for cytokine 1, anti-FITC antibody green fluorescence conjugate, capture antibody for cytokine 2, biotinylated detection antibody for cytokine 2, streptavidin-phycoerythrin conjugate, BSA. Kit can be sold with or without plates. |
| Tested Applications | Fluorospot                                                                                                                                                                                                                                                                                                                                     |
| Brand               | Diaclone                                                                                                                                                                                                                                                                                                                                       |
| Storage             | Store at +4°C                                                                                                                                                                                                                                                                                                                                  |
| Cross Reactivity    | No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10                                                                                                                                                                                   |
| Specificity         | Recognizes natural human IL-10 and IFN-g                                                                                                                                                                                                                                                                                                       |