

Product Datasheet

Human IFNγ / IL-10 Dual Fluorospot Set 874.012.015

Product Name	Human IFNγ / IL-10 Dual Fluorospot Set
Catalog Number	874.012.015
Description	Diaclone Dual Fluorospot Sets include capture antibody for cytokine 1, FITC-conjugated detection antibody for cytokine 1, anti-FITC antibody green fluorescence conjugate, capture antibody for cytokine 2, biotinylated detection antibody for cytokine 2, streptavidin-phycoerythrin conjugate, BSA. Kit can be sold with or without plates.
Tested Applications	Fluorospot
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-γ. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10
Specificity	Recognizes natural human IL-10 and IFN-γ