

## Product Datasheet

### Human IFN $\gamma$ / IL-10 Dual Fluorospot Set 874.012.020

<b>Product Name</b>	Human IFN $\gamma$ / IL-10 Dual Fluorospot Set
<b>Catalog Number</b>	874.012.020
<b>Description</b>	Diaclone Dual Fluorospot Sets include capture antibody for cytokine 1, FITC-conjugated detection antibody for cytokine 1, anti-FITC antibody green fluorescence conjugate, capture antibody for cytokine 2, biotinylated detection antibody for cytokine 2, streptavidin-phycoerythrin conjugate, BSA. Kit can be sold with or without plates.
<b>Tested Applications</b>	Fluorospot
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10
<b>Specificity</b>	Recognizes natural human IL-10 and IFN-g