

Product Datasheet

Anti-h Copeptin 4802 SPTN-5 100639

Product Name	Anti-h Copeptin 4802 SPTN-5
Catalog Number	100639
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Epitope is located within the sequence ATQLDGPAGALLRLVQ in the N-terminal part of copeptin.
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
IEF Profile	6.4-6.8
Cross Reactivity	N/D
Specificity	Antibody recognizes human copeptin
Shelf Life	Unspecified