

Product Datasheet

Anti-h Copeptin 4804 SPTN-5  
100649

|                     |  |
|---------------------|--|
| Product Name        | Anti-h Copeptin 4804 SPTN-5  |
| Catalog Number      | 100649   |
| Description         | Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components. |
| Tested Applications | FIA  |
| Brand               | Medix Biochemica   |
| Form/Appearance     | Liquid, may turn slightly opaque during storage  |
| Concentration       | 5.0 mg/ml (+/- 10 %)   |
| Storage             | +2-8°C   |
| Isotype             | IgG1   |
| Clonality           | Monoclonal   |
| Epitope             | Epitope is located within the sequence QLAGAPEPFEPAPDAY in the C-terminal part of copeptin.        |
| Purity              | ≥ 95 %   |
| Buffer              | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative                               |
| IEF Profile         | 7.2-7.8  |
| Cross Reactivity    | N/D  |
| Specificity         | Antibody recognizes human copeptin   |
| Shelf Life          | Unspecified  |