

Product Datasheet

Anti-Human IL-6 Capture Antibody, Azide Free Clone B-E8 879.030.001

Product Name	Anti-Human IL-6 Capture Antibody, Azide Free Clone B-E8
Catalog Number	879.030.001
Description	This antibody can be used as Capture in a human sandwich Immunoassay. . Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, ELISpot
Immunogen	Recombinant human IL-6
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,5-2,5µg/ml ELISpot: 5-10µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises both natural and recombinant human IL-6
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells