

Product Datasheet

Anti-Human IL-7 Capture Antibody, Azide Free Clone B-N18 879.040.001

Product Name	Anti-Human IL-7 Capture Antibody, Azide Free Clone B-N18
Catalog Number	879.040.001
Description	This Antibody can be used as Capture in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human IL-7
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,5-2,5µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant IL-7
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells