

Product Datasheet

Anti-Human IL-10 Capture Antibody, Azide Free Clone B-S10 879.100.001

Product Name	Anti-Human IL-10 Capture Antibody, Azide Free Clone B-S10
Catalog Number	879.100.001
Description	This Antibody can be used as Capture in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, ELISpot
Immunogen	Recombinant human IL-10
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 1-5µg/ml ELISpot: 5-10µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human IL-10
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells