

Product Datasheet

Anti-Human CD25 Detection Antibody, Biotin Conjugated Clone B-F2 879.500.002

Product Name	Anti-Human CD25 Detection Antibody, Biotin Conjugated Clone B-F2
Catalog Number	879.500.002
Description	This biotinylated antibody can be used as Detection in a human sandwich Immunoassay. . Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	PHA activated T cells
Alternative Names	IL-2R
Brand	Diaclone
Form/Appearance	Liquid
Concentration	0.1 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,05-0,5µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	PHA activated T cells
Buffer	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Specificity	Recognises the Interleukin 2 Receptor a chain (IL-2 Ra), a 55 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells