

Product Datasheet

Anti-Human CD126 Capture Antibody, Azide Free Clone B-N12 879.510.001

Product Name	Anti-Human CD126 Capture Antibody, Azide Free Clone B-N12
Catalog Number	879.510.001
Description	This Antibody can be used as Capture in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	gp80 transfected CHO cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,3-1µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Interleukin-6 Ra subunit, IL-6 Ra, IL-6 R, a 80 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells