

Product Datasheet

Anti-Human CD130 Detection Antibody, Biotin Conjugated Clone B-T2 879.520.002

Product Name	Anti-Human CD130 Detection Antibody, Biotin Conjugated Clone B-T2
Catalog Number	879.520.002
Description	This biotinylated Antibody can be used as Detection in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Natural soluble gp130
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. With 1% BSA and 0,09% Sodium Azide
Storage	Stable at +2-8°C for 12 months.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,05-0,5µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Gp130, common subunit for IL-6, IL-11, OSM, LIF, CNTF, CT-1 receptors, a 130-140 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells