

Product Datasheet

Anti-Human CD62E Detection Antibody, Biotin Conjugated Clone B-S3 879.530.002

Product Name	Anti-Human CD62E Detection Antibody, Biotin Conjugated Clone B-S3
Catalog Number	879.530.002
Description	This biotinylated Antibody can be used as Detection in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Activated human endothelial cells
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. With 1% BSA and 0,09% Sodium Azide
Storage	Stable at +2-8°C for 12 months.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,05-0,5µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Endothelial Leucocyte Adhesion Molecule-1 (ELAM-1), E-selectin, a 97 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells