

Product Datasheet

Anti-Human p53 Capture Antibody, Azide Free Clone B-F34 879.630.001

Product Name	Anti-Human p53 Capture Antibody, Azide Free Clone B-F34
Catalog Number	879.630.001
Description	This Antibody can be used as Capture in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant p53 protein
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Store at +2-8°C for 12 months. DO NOT FREEZE.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 1-5µg/ml
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognises wild type and mutant p53, the epitope is within residues 18-30 of human p53
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells