

Product Datasheet

Anti-Human p53 Detection Antibody, Biotin Conjugated Clone B-P3 879.630.002

Product Name	Anti-Human p53 Detection Antibody, Biotin Conjugated Clone B-P3
Catalog Number	879.630.002
Description	This biotinylated Antibody can be used as Detection in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant p53 protein
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. With 1% BSA and 0,09% Sodium Azide
Storage	Stable at +2-8°C for 12 months.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,05-0,5µg/ml
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises wild type and mutant p53, the epitope is within residues 18-30 of human p53
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells