

Product Datasheet

Anti-Human CD141 Detection Antibody, Biotin Conjugated Clone B-S24 879.720.002

Product Name	Anti-Human CD141 Detection Antibody, Biotin Conjugated Clone B-S24
Catalog Number	879.720.002
Description	This biotinylated Antibody can be used as Detection in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Eahy 926 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. With 1% BSA and 0,09% Sodium Azide
Storage	Stable at +2-8°C for 12 months.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,05-0.5µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Recombinant antibody from the original clone
Specificity	Recognises thrombomodulin (TM), a 105 kDa protein