

Product Datasheet

Anti-Human IL-23 Capture Antibody, Azide Free Clone B-Z23 879.920.001

Product Name	Anti-Human IL-23 Capture Antibody, Azide Free Clone B-Z23
Catalog Number	879.920.001
Description	This Antibody can be used as Capture in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant IL-23 (sub-unit p19)
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 2,5-5µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant IL-23, p19 sub-unit
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells