

Product Datasheet

Human IL-4 ELISA Kit 950.020

Product Name	Human IL-4 ELISA Kit
Catalog Number	950.020
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	1.1-35 pg/ml
Sensitivity	0.7 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-4
Sample Types	Serum, Plasma, Cell culture supernatant