

Product Datasheet

Anti-h ST2 10204 SPTN-5 100683

Product Name	Anti-h ST2 10204 SPTN-5
Catalog Number	100683
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	CLIA, FIA
Alternative Names	Soluble Interleukin 1 Receptor-like 1
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Recognizes free sST2
Purity	≥ 95 %
Affinity constant	N/D
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	5.6-6.2
Cross Reactivity	Not determined
Specificity	Antibody recognizes human soluble ST2
Shelf Life	36 months