Medix Biochemica

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 Fl-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

Human IL-6 ELISA Kit 950.030

Product Name Human IL-6 ELISA Kit

Catalog Number

950.030

Description

Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed

procedure.

Range 6.25-200 pg/ml

Sensitivity 2 pg/ml

Brand Diaclone

Storage Store at +4°C

Cross Reactivity No cross reactivity with other human cytokines. No interference

with soluble receptor gp80.

Specificity Recognizes both natural and recombinant human IL-6

Sample Types

Serum, Plasma, Cell culture supernatant