

Product Datasheet

Human IL-6 HS ELISA Kit 950.035

Product Name	Human IL-6 HS ELISA Kit
Catalog Number	950.035
Description	Diaclone High Sensitivity ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	1.56-50pg/ml
Sensitivity	0.81 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. No interference with soluble receptor gp80.
Specificity	Recognizes both natural and recombinant human IL-6
Sample Types	Serum, Plasma, Cell culture supernatant