

Product Datasheet

Human TNF alpha ELISA Kit 950.090

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Human TNF alpha ELISA Kit |
| Catalog Number | 950.090 |
| Description | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range | 25-800 pg/ml |
| Sensitivity | 8 pg/ml |
| Tested Applications | ELISA |
| Alternative Names | TNFa |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with soluble TNFR I and TNFR II. |
| Specificity | Recognizes both natural and recombinant human TNF-a |
| Sample Types | Serum, Plasma, Cell culture supernatant |