

# Product Datasheet

## Human TNF alpha ELISA Kit 950.090

Product Name	Human TNF alpha ELISA Kit
Catalog Number	950.090
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	25-800 pg/ml
Sensitivity	8 pg/ml
Tested Applications	ELISA
Alternative Names	TNFa
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with soluble TNFR I and TNFR II.
Specificity	Recognizes both natural and recombinant human TNF-a
Sample Types	Serum, Plasma, Cell culture supernatant