

## Product Datasheet

### Human TNF alpha ELISA Kit 950.090

<b>Product Name</b>	Human TNF alpha ELISA Kit
<b>Catalog Number</b>	950.090
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	25-800 pg/ml
<b>Sensitivity</b>	8 pg/ml
<b>Alternative Names</b>	TNF-alpha, TNF-a
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with soluble TNFR I and TNFR II.
<b>Specificity</b>	Recognizes both natural and recombinant human TNF-a
<b>Sample Types</b>	Serum, Plasma, Cell culture supernatant